

NOBAMED Paul Danz AG

NOBAFIX® Individually packed in foil

REF 410380

Product Description, Intended use, Application

This fixation bandage with a width of 8 cm and a stretched length of 4 m is a rolled, lengthwise-elastic dressing retention bandage. Due to its weave, the bandage has a fine, smooth and tight structure. The elastic bandage is used to fix wound pads, especially over joints and body curves. As a damp bandage, NOBAFIX® is also used to shape and/or fix plaster splints and cast dressings. NOBAFIX® is for single use.

Composition

Polvester

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PP foil

Secondary packaging: folding box made of

cellulose

Tertiary packaging: carton made of

cellulose

Symbols used in labelling

Explanations at www.nobamed.com







Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

Storage and Transport

Dry and dustfree, protect from sunlight.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: 03.05.2021 [REV 7], replaces: 11.05.2020 [REV 6]