

310250 Mefix

Self-adhesive fabric

Product details

Size :	2,5cm x 10m
Descriptive feature :	Self-adhesive fixation
Sterile :	Non-sterile
Brand :	Mefix

Images



Delivered items

310250-30

Sales released in: Algeria, Austria, Bahrain, Belgium, Brazil, Bulgaria, Canada, Chile, China, Croatia, Cyprus, Czechia, Denmark, Estonia, Faroe Islands, Finland, France, Germany, Greece, Hong Kong, Hungary, Iceland, Iran (Islamic Republic of), Iraq, Ireland, Italy, Japan, Kazakhstan, Korea (the Republic of), Kuwait, Latvia, Lebanon, Lithuania, Luxembourg, Malaysia, Malta, Mauritius, Netherlands, Norway, Oman, Pakistan, Poland, Portugal, Qatar, Romania, Saudi Arabia, Serbia, Singapore, Slovakia, Slovenia, South Africa, Spain, Sweden, Switzerland, Taiwan (Province of China), Thailand, Ukraine, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland, Viet Nam

Country of origin: Thailand

Shelf life: 3 years

Sterilization method: Non-sterile

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

Packing information: First packaging layer is a cardboard dispenser box. Second layer is a corrugated board

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2023 Mölnlycke Health Care AB. All rights reserved.



transport box.

Is suitable for Tray: Yes

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430082282			
Consumer pack	1	7323190203246			
Dispenser box	1	7323190195107	31x74x72 mm		
Transport box	40	7323190195091	195x263x157 mm	8.1 dm3	2 / 1.5 kg
Pallet	7920	7323190195084	800x1200x1871 mm	1796.2 dm3	421.0 / 301.0 kg

Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Product Composition Film Products and Fixation Products

Product Component	Composition
Backing material	Hydrophobic polyester nonwoven
Adhesive	Polyacrylate
Protective release liner	Siliconised paper

Product Performance Fixation Products

Characteristics	Test Method	Internal Test Method	Unit	Requirement	Product Performance
Adhesion	ASTM D3330/D3330M-04, method F	T-1082	N/25 mm	Not specified	2.97
Moisture Vapour Transmission Rate (MVTR) of a wound dressing when in contact with water vapour	EN 13726-2 part 3:2	T-1070	g/m ² /24 h	Not specified	3448

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2023 Mölnlycke Health Care AB. All rights reserved.



Technical

Dimension

Dimension text	Dimension value
Product	1 in x 11 yd
Product	2.5 cm x 10 m

Classifications

Regulation type(s)	MDR Class I ns	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :	1			
Conformity Annexes :	IV			
Measuring Function :	No			
Notified body medical devices/PPE :	MPA	Local Authority Thailand	Local Authority EU Non EEC	Local Authority Thailand
510(k) clearance number :		NA		
Regulatory Released :	Austria, Belgium, Bulgaria, Croatia, Cyprus, Czechia, Denmark, Estonia, Faroe Islands, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland	Thailand	Serbia	Thailand

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2023 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority MEA	Local Authority South Korea	Local Authority Taiwan	Local Authority United Arab Emirates
510(k) clearance number :				
Regulatory Released :	Algeria, Bahrain, Lebanon, Mauritius, Oman, South Africa, Ukraine	Korea (the Republic of)	Taiwan (Province of China)	United Arab Emirates

Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority China	Local Authority MEA	Local Authority Pakistan	Local Authority Americas
510(k) clearance number :		NA		
Regulatory Released :	China	Qatar	Pakistan	Chile

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2023 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Saudi Authority	Local Authority Kuwait	Local Authority APAC	Local Authority MEA
510(k) clearance number :				NA
Regulatory Released :	Saudi Arabia	Kuwait	Hong Kong, Malaysia, Singapore	Algeria, Bahrain, Lebanon, Mauritius, Oman, South Africa, Ukraine

Classifications

Regulation type(s)	Locally Regulated	Class I	Class I	Non Medical Device
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority CIS	Health Canada	Anvisa (Brazil)	PMDA (Japan)
510(k) clearance number :	NA			
Regulatory Released :	Kazakhstan	Canada	Brazil	Japan

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN ISO 20417, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 10993-5, EN ISO 10993-10, EN ISO 10993-18, ISO 14001

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2023 Mölnlycke Health Care AB. All rights reserved.



Removable label

No

GMDN Code (Global Medical Device Nomenclature)

16866 Basic adhesive tape

UNSPSC

42311552 Medical and surgical adherent tapes for general use

Commodity Code

3005100000 Wadding, Gauze, Dressings, drapes singlepacked - adhesive articles, Sets mainly consisting thereof

CE Responsibility

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

733243000000000012JF

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2023 Mölnlycke Health Care AB. All rights reserved.

